




3D custom-made implants



Sebbin
PARIS

3D CUSTOM-MADE IMPLANTS

For Correction of Pectus Excavatum and Poland Syndrome



An Innovative Partnership



Specialist in designing 3D custom-made implants, after 10 years of collaboration with the University Hospital in Toulouse

Scientific support of Pr. Chavoin




Experts in the conception and manufacturing of silicone gel implants

Working in partnership with surgeons for patient's care in plastic and aesthetic surgery


**A 100 % Personalized Solution
Perfectly suitable for deformities such as
Pectus Excavatum and Poland Syndrome**

Thoracic Deformities



Pathologies - Pectus Excavatum

- Congenital thoracic deformity characterised by a medial or lateral depression of the sternum
- It is the most common congenital thoracic malformation
- Very rare functional troubles, with major psychological impact



1 to 2%

of people affected

Only 15%

of people affected by the Pectus Excavatum undergo surgery

Occuring in **1 over 300**
births

Pathologies - Pectus Excavatum


3 TYPES ACCORDING TO CHIN

Type 1
Symmetrical, deep and
focuses on the sternum



Pectus Excavatum Type 1 in a Woman

Type 2
Symmetrical, shallow,
extends to the pectoral regions



Pectus Excavatum Type 2 in a Man

Type 3
Asymmetrical and extends to the pectoral
regions deviation is mostly on right side



Pectus Type 3 in a Man

OTHER TYPES



Ravitch Revision in a Man

Pectus Arcuatum
Mixed Pectus
Revisions



Pectus Arcuatum in a Woman

Pathologies - Poland Syndrome

- Relatively rare malformation
- Complete or partial absence of the pectoralis muscle
 - Frequent occurrence : Depression under the clavicle and mammary asymmetry in women
- Can be associated with a homolateral hand deformity


Occuring in **1 over 30 000**
births



Pathologies - Poland Syndrome

3 TYPES

Type 1
Correct subcutaneous plane



Poland type 1 in a Man

Type 2
More atrophic subcutaneous plane



Poland type 2 in a Woman

Type 3
Atrophy of subcutaneous plane with
thoracic deformity





Poland type 3 in a Woman




Poland type 2 in a Man

Pathologies - Mammary Deformities


- The Pectus Excavatum and the Poland Syndrome affect the mammary shape in women (convergence, asymmetry, divergence...)
- The thoracic deformity should be treated primarily
- Mammary implants can be placed in a secondary operation if necessary (minimum 6 months later)



Pectus Type 1
Strong Convergence



Pectus Type 2
Slight Convergence



Pectus Type 3
Asymmetry


Custom-Made Solution



Precise Technology - Conception



3D Reconstruction of patient's body from a CT scan



Virtual aided design of the implant using the patients anatomy. Taking into account anterior anatomical plane and posterior surgical plane of the thorax




Control and optimisation of the result with the surgeon

Precise Technology - Manufacturing



- Manufacturing of the implant's prototype from the 3D design
- Manufacturing of a mould in which, the silicone elastomer will be injected






A Simple and Minimally Invasive Surgery

The Surgical Technique (*Pectus Excavatum*)

- Pre-operative marking, utilizing the prototype and its anatomical markers
- Incision, vertical median - 7 cm length
- Dissection in retro-muscular (pectorals) and retro-aponevrotic (abdominals)
- Preparation of the cavity to fit the implant
- Insertion of the flexible implant in deep position
- Suture in 3 planes



The surgery time is generally **1 hour**

1 to 2 days of hospital stay

A Simple and Minimally Invasive Surgery



Post Surgery

- Moderate pain
Treatable with simple painkillers
- Low post-operative care
Suture, no drain
Dressing maintained 8 days
- Compression bra
Duration 1 month (24/7)
- Follow-up consultation every 8 days
With drainage puncture of seroma if necessary
- 15 Days off-work
Progressive return to sport, not before 3 months



Contention bra



Puncture of seroma

Rapid Recovery



Lifelong Implants

Medical Grade Silicone Elastomer


- Semi-rigid rubber
No risk of rupture or retraction
No risk of rejection, encapsulation and stabilisation of the implant
- Smooth Surface
No cases of lymphoma reported



Lifelong implants, no need of replacement.

Immediate Aesthetic Results

Pectus Excavatum



Pectus Type 1 in a Woman



Pectus Type 1 in a Man



Pectus Arcuatum in a Man



Pectus Type 3 in a Woman



Nuss Revision in a Woman




Ravitch Revision in a Man

Immediate Aesthetic Results

Poland Syndrome



Poland type 1 in a Man



Poland type 3 in a Woman



Poland type 2 in a Man



Poland type 2 in a Woman

Custom-Made Implants

A good alternative to traditional orthopedic techniques

In the absence of functional troubles

In complex cases of Asymmetrical Pectus or Pectus Arcuatum

For secondary surgeries after failure of Nuss or Ravitch

For all patients from 15 to 65 years old

Functional tests can be performed to validate the choice of the technique



A Simple Process



1 CT SCAN

The patient's entire chest is scanned

- Supine with arms along the body
- 1 to 1.2 mm thick DICOM standard-format cuts
- Medical tomography with a CT scan or mediastinal filter to achieve a good contrast

2 PATIENT SCAN AND PRESCRIPTION SENT

3 QUOTE DRAWN UP



4 QUOTE APPROVED AND ORDER CONFIRMED

5 3D IMPLANT DESIGNED AND APPROVED BY THE SURGEON

6 IMPLANT MANUFACTURED AND STERILISED




7 DELIVERY

8 to 12

weeks from the validation of the implant design by the surgeon

Literature Available

- 1- Chavoin J.P., et al. Correction of Pectus Excavatum by Custom-Made Silicone Implants: Contribution of Computer-Aided Design Reconstruction. A 20-Year Experience and 401 Cases. *Plast Reconstr Surg.* 2016.
- 2- Chavoin J.P., et al. Correcting Poland Syndrome with a Custom-Made Silicone Implant: Contribution of Three-Dimensional Computer-Aided Design Reconstruction. *Plast Reconstr Surg.* 2018.
- 3- Chavoin, J.P., (Ed.) *Pectus Excavatum and Poland Surgery. Custom-Made Silicone Implants by Computer Aided Design.* Springer. 2019.
- 4- Chavoin J.P., et al. The Role of Computer-Aided Design Implant Insertion in Revision Pectus Surgery. *The Annals of Thoracic Surgery* Volume 112, Issue 5, November 2021, Pages e387-e390.



A Large Surgeon Network

Around 200 referral surgeons worldwide

- Thoracic, Plastic, Pédiatric surgeons
- Trained to the technique of custom-made implants
- Personal consultation with Pr Chavoin
- Visible on AnatomikModelings website (12k visitors a month) and are able to be directly contacted by patients





3D custom-made implants



Sebbin
PARIS